

# BRI

**Mind and Body Technology**

*Foster Self Knowledge.  
Improve Lives.*



*Biofeedback/ Neurofeedback Equipment, Supplies  
BCIA Certification Training and Consulting*

**Phone:** (877) 669-6463 / (914) 762-4646

• **Fax:** (914) 762 2281

• **Email:** [info@biofeedbackinternational.com](mailto:info@biofeedbackinternational.com)

• **Visit our Website:** [www.biofeedbackinternational.com](http://www.biofeedbackinternational.com)

109 Croton Ave  
2nd Floor, Suite #240  
Ossining, New York 10562



*Foster Self Knowledge.  
Improve Lives.*

## **MISSION**

**Foster self knowledge and improve lives through innovation, continuous improvement of products, services and employee productivity, and relentless efforts to deliver highest values at lowest prices.**

## **WHO ARE WE?**

**Biofeedback Resources International Corporation (BRI) was founded by R. Adam Crane in order to meet the needs of the growing biofeedback community and the emerging public market worldwide. Since 1971, Mr. Crane and the teams he has assembled have brought many innovations and services to the biofeedback community. After 22 years of working together Mr. Crane sold the business to Harry L. Campbell.**

**Our BRI core team includes individuals who have outstanding reputations in the biofeedback community for their skill, integrity and dedication.**

**We want to be your biofeedback company, and assert that we can assist you in succeeding as well or better than any company in our field, because we are leaders in many critical areas of Biofeedback and Neurofeedback. A visit to our website will give you insight into our strengths and innovations.**

## **TRAINING (See BRI Training Catalog)**

**Our training Division, Health Training Seminars, offers at many comprehensive programs every year which are accredited by the Biofeedback Certification Institute of America and the American Psychological Association as well as specialized biofeedback training, online courses, mentoring, coaching and consulting (in person and online). Our programs offer some of the most comprehensive and cost effective biofeedback and neurofeedback training in the world.**

---

## Table of Contents

---

Computerized Biofeedback/EEG Neurofeedback Systems	4
Stand Alone Instruments	8
Microcurrent Pain Management	13
Computer Software-General /Application	14
Supplies and Electrodes	17
Disposable Sensors	18
J&J Supplies	19
Continance Supplies	20
Thought Technology Sensors	21
Audio Cassettes, CDs, DVDs, and VHS Tapes	22
Books	25
Computer Peripherals	29
Computer Furniture	29
Customer Support, Consulting, Training	29
Terms & Conditions	29

**IMPORTANT NOTE: PLEASE VISIT OUR WEBSITE:**

**WWW.BIOFEEDBACKINTERNATIONAL.COM FOR THE LATEST UPDATES ON TRAINING  
PROGRAMS AND EQUIPMENT.**

**For special offers and discounts, please call us toll free at 877-669-6463**

---

Thought Technology Systems

**PROCOMP INFINITI:**

The ProComp Infiniti is an 8 channel, multi-modality encoder that has all the power and flexibility you need for real-time, computerized biofeedback and data acquisition in any clinical setting. Includes three sets of software: Biograph Infiniti, Choice of Physiology or EEG Suite.




Modalities include: EMG, EEG, Temp, Skin Conductance, Heart Rate/BVP, Respiration, and EKG



ProComp Infiniti comes complete with:

- 14 bit resolution, eight channel
- ProComp Infiniti encoder unit
- TT-USB interface unit
- Fiber-optic cables (1' and 15')
- Four alkaline "AA" batteries
- Sleek fabric storage and carrying case

See Supplies and Electrodes for available sensors



Price: \$3995.00  
Order #: T7500M  
Sensors additional

**PROCOMP 2**

The ProComp 2 is a compact yet powerful 2 channel device that allows clinicians to set up a second clinical system—or to empower their clients by offering a take-home device that is convenient to wear on a headband or a shirt collar. The ProComp 2 contains a built-in EEG sensor (simply connect an extender cable for EEG monitoring and biofeedback), and it can use any of the two ProComp Infiniti sensors. The ProComp 2 systems contains all the peripherals to easily connect it to a desktop or laptop IBM-compatible PC.

The ProComp 2 with the powerful Physiology or EEG Suite Software

Price: \$1850.00  
Order #T7400M

ProComp 2 comes complete with:

- 13 bit resolution, 2 channel ProComp 2 encoder unit
- Fiber-optic cables ('1 and 15')
- Two alkaline "AA" batteries
- Carrying case



ProComp5 Infiniti™ is housed in an ergonomically-designed case and requires only a USB port to connect to any IBM compatible PC. ProComp5 Infiniti™ has the identical inputs as ProComp Infiniti, but omits the last 3 channels. The first two sensor channels provide ultimate signal fidelity (2048 samples per second) for viewing raw EEG, SEMG, EKG and HR/BVP signals, while the remaining 3 channels sample data at 256 samples/second for slower signals such as respiration, temperature, force, etc. Software includes Biograph Infiniti and either the EEG Suite or the Physiology Suite.

Price: \$3,150.00  
Order #: T77525  
Sensors additional

## Computerized Systems

### Thought Technology Systems

#### **FLEXCOMP INFINITI**

Our top-of-the-line encoder is the ideal data acquisition and physiological monitoring device for power users. It offers 10 high speed channels (2048 samples/seconds) with 14 bits of resolution (1 part in 16364) and acquire data from any Thought Technology sensor.

Researchers, SEMG specialists and ergonomists will appreciate it's ability to monitor and record SEMG and displaying it as raw , median frequency, RMS, peak-to-peak, 2D and 3D frequency and power spectrum. Up to four encoders (each using 10 channels) can be connected simultaneously to produce the largest possible "window" into the subject's physiology.

Flexcomp Infiniti offers, internal, user-activated calibration to ensure that you can always obtain the highest quality signal, without the costly downtime associated with factory re-calibration.

Flexcomp Infiniti comes complete with: Price: \$6,495.00

- 14 bit resolution, 10 channel Flexcomp Infiniti encoder unit Order #T7555M

- TT-USB interface Unit
- Fiber-optic cable ('1 and 15')
- Four alkaline "AA" batteries
- Sleek fabric storage and carrying case

#### J & J Systems

Price \$1,995  
Order# I-330 C2+

#### **I-330 C2+ 12 Channel**

The popular I-330-C2+, now has advanced Windows software and a fast USB interface! It's 4 differential channels for EMG, ECG, EEG, and 8 channels for temperature, skin resistance, and respiration, provide power, flexibility and functionality for a multitude of customized uses. J&J Engineering's new Windows software displays signals, provides feedback, collects data, prints reports, and exports database compatible files.

Features include automatic testing of electrode impedance for connection quality, easy hookup with reusable gel free sensors, and fast channel sampling at 1024 SPS for high resolution spectral displays and precision filtering. Connects via USB to any 1 GHz (or faster) laptop or desktop running WindowsXP or older Windows versions with current updates.

Features:

\* 2-person monitoring for ECG, HRV, Respiration, Skin Resistance or Conductance, and Temperature

\* 4 EMG channels with spectral and raw signal displays, 4 Skin Resistance and Conductance, 4 Temperature

\* ECG, 2 EMG, 2 Respiration, 2 Skin Resistance, 2 Temperature

\* 1, 2, or 4 EEG channels, Respiration, SR, Temperature

**For special offers and discounts, please call us at 877-669-6463**



#### **C-2+ OPTIONS**

SE-22A Ear Clip, each  
\$20

**MC-6SY:** C2+ Temperature and skin resistance cable  
**\$ 75.00** #MC-6SY

**MC-5SGW40:** C2+ Electrode cable for EMG/EEG/  
ECG Includes TP-60 and TP-40 Snap Leads  
**\$ 80.00** #MC-5SWG


**MC-5SGW18:** C2+ Electrode cable for EMG/EEG/  
ECG Includes TP-60 and TP-18 Snap Leads  
**\$ 75.00** #MC-5SWG

**TP-60** C2+12 Electrode cable with three 1.5mm jacks  
for TP leads or standard EEG electrodes  
**\$ 60.00** #TP-60

**MC-3MY:** C2+ Respiration Sensor  
**\$ 150.00** #MC-3MY

**MC-5D:** C2+ Gel Free Wrist Electrodes Optional  
Reusable Heart Rate Sensor  
**\$ 85.00** #MC-5D

## Computerized Systems



**SE-35:** C2+ EDR/SR Velctro Finger Electrodes Set of 2

**\$ 12.00 #SE-35**

**SE-13:** C2+ Replacement Electrodes for EEG Ear-clip and ECG Wrist Electrodes, Pack of 6

**\$ 10.00 #SE-13**

**SE-12:** C2+ Replacement Electrodes for SE-20 EEG Cups, Pack of 100

**\$ 15.00 #SE-12**

**MC-EXT:** C2+ 12 Cable Extender Adds 6 ft to MC type cables

**\$ 85.00 #MC-EXT**

**SE-55:** C2+12 Pre Gelled Disposable Electrodes 1 3/8th" Pack of 100

**\$ 38.00 #SE-55**

**SK-1:** EEG Supplies Kit Includes SE-29 Electrodes and SE-22A earclips

**\$ 60.00 #SK-1**

**SE-20:** C2+ Disposable EEG Electrodes Pack of 100 cups and pellets

**\$ 30.00 #SE-20**




The CapnoTrainer is a USB-based, software-operated, educational capnograph for learning about breathing habits and their effects. The instrument is for teaching and learning new breathing behavior patterns. The instrument is not intended for diagnosis and/or treatment.

**Price \$3,300 #CapnoTrainer-M**

CapnoPlus ECG, EMG, HRV Add-on Price \$850 #CapnoPlus

### emWave Desktop

is proving to be a huge asset to those who are ready to add a new level of emotional mastery to their skills.



The EmWave PC is a simple, easy to use software program that displays your heart rhythms and shows you how stress may be affecting you. In this program, you'll learn two tools, Freeze-Frame® and Heart Lock-In® to help make your life easier in a busy and

What would you like to improve? Your work performance; your school work; your golf game; your creative efforts; your health; your fitness; your work; home or social relationships; your peace of mind?

All of these things are affected by the overwhelm, the pace and intensity of life we know as stress and your ability to deal with it. Let the emWave PC help you break through stress that's preventing you from feeling good, being in control and achieving your goals.

Price \$295  
Order# EMWAVEPC

**For special offers and discounts, please call us at 877-669-6463**

## Thought Technology Systems

### **PROCOMP 2**

#### **Biofeedback Starter System - Stress Control Suite**


##### Quick Overview

Our starter biofeedback systems are ideal for clinicians who want to begin exploring physiological biofeedback without overwhelming their budget. This medical grade system provide you with out-of-the-box functionality and the possibility of expanding your system when you are ready to go beyond the basics.

##### Details

Introduce biofeedback into your practice with this compact and easy to use two-channel device. The system, which includes skin conductance and temperature sensors, comes with powerful multimedia biofeedback software and the Stress Control Suite.

Price: \$2,355  
Order #: T7405



### **PROCOMP 2**

#### **HRV Starter System**

##### Quick Overview

Our starter HRV system is ideal for clinicians who want to begin exploring heart rate variability biofeedback without overwhelming their budget. This medical grade system provides out-of-the-box functionality with a set of simple tools for assessing and training HRV.

##### Details

Introduce heart rate variability biofeedback into your practice with this compact and easy to use two-channel device which comes with a finger pulse and a respiration sensor, powerful multimedia biofeedback software and the HRV Suite.

Perfect for psychologists and other therapists who want to discover the benefits of heart rate variability biofeedback with simple, yet powerful tools. With the HRV suite, you can perform a basic HRV assessment, run biofeedback sessions and teach your client deep relaxation skills with paced breathing exercises.

Price: \$2,370  
Order #: T7406



#### **MyoTrac 3 w/ Rehab Suite - T9920**

The MyoTrac3 with rehab suite is a highly sensitive, portable, computerized, dual-channel sEMG interface. The Windows based continence protocol software allows clinicians assessment flexibility and the ability to customize applications.

Price: \$1,995



## Computerized Systems



**BRAINMASTER Atlantis I (4x4):** Atlantis I (4+4) goes beyond the state of the art to define a new approach to high quality and effective EEG and peripheral feedback. It includes an innovative Triamp (tm) bioamplifier that records and reports the EEG signal, the DC and slow cortical potentials, and the impedances of the sensors as well, all continuously, and in real time. In addition to 4 (24-bit) EEG channels, there are 4 (12-bit) Aux channels, also using Triamp (tm) inputs, ready to provide simultaneous AC, DC, and impedance recording for skin, muscle, heart, temperature, and other add-on peripheral modalities expected in the near future.

Price: \$2,695.00  
Order# ATLN-1



Price: \$1,995.00  
Order# ATLN-2



**BRAINMASTER Atlantis II (2x2):**

ATLANTIS II is a 2x2 design (with two channels of EEG and two channels of AUX signals) ideal for additional biofeedback modalities. Includes continuous real-time impedance checking and recording; Total Immersion with photic, vibratory and auditory feedback. Evolving 24 bit USB technology.

ATLANTIS I and II: Biofeedback with Impedance checking; USB; Direct Current (DC) and Slow Cortical Potential (SCP) measurement, using Triamp amplifier design.

**The Discovery 24E** is a 24-channel EEG and DC amplifier that is an evolution of the technology used in the 2E and 2EA (Atlantis) series of EEG biofeedback devices. It is currently available for field testing. It is fully compatible with QEEG systems such as NeuroGuide , Loreta, SKIL , and WinEEG, and incorporates BrainMaster training capabilities including amplitude and connectivity training, including Live Z-Score training. \$5,800 #Discovery-24



**Blood Pressure: DynaPulse Pathway:** the Pathway System uses software patterned after the DynaPulse 200M series. Though this software has a newly designed algorithm that produces blood pressure measurements that correlate to readings obtained manually in a physician's office. The Pathway will track unlimited number of patients and their data while producing a choice of trend graphs to report the blood pressure measurements. It is an auto-pump for easy cuff inflation and is AC/DC powered.

Price: \$995  
Order# DP-200M



## Stand Alone Instruments

---

### *Electromyographs—EMG Muscle Biofeedback*

#### MYOTRAC

Myotracs packs extreme sensitivity into a light, compact unit that works for you and your patient in the office or home.

##### Features:

- The highly-sensitive Myoscan active sensor needs no skin preparation
- Two modes of biofeedback –LED and 5 way tone - provides flexibility
- 5 way tone feedback includes proportional, above and below threshold and 4 second delayed alarm.
- 13-segment LED bargraph offers expanded range for precise visual feedback
- Wide and narrow bandpass filters provide precise monitoring of normal and fatiguing muscles and eliminate heartbeat artifact
- 0 to 20 volts DC output allows to connect to third party equipment and computerized monitoring
- Low current consumption provides long battery life



Price \$499  
Order #T4000P

MyoTrac Infiniti systems allow you to measure high-resolution surface electromyography (SEMG) and perform electrical stimulation in a large number of configurations. The MyoTrac Infiniti stands apart from the competition by also performing SEMG-triggered stimulation. SEMG is used to trigger stimulation, to help patients improve volitional activation of target muscle groups. This combination of three modalities enhances treatment possibilities and eliminates the requirements for multiple devices.



Price: \$2,500  
Order#: T9850US  
Clinical w/Rehab SW

Price: \$1,600  
Order#: T9800US  
Stand-Alone

For special offers and discounts, please call us at 877-669-6463

---

## Stand Alone Instruments

### Electromyographs—EMG Muscle Biofeedback

#### U-CONTROL EMG

Enables patients to self-train with prescribed exercises, at home to, to re-educate pelvic muscles for better bladder and bowel functions.



UCONTROL has an adjustable threshold setting to enable health professionals and patients to revise and report on their training programs as they progress.

Patients use a comfortable, reusable vaginal or rectal sensor in combination with simple biofeedback exercises to improve bladder and bowel functions by:

- Holding a steady LED readout or tone
- Intensifying muscle contractions above the threshold level
- Repeating a pre-set number of exercises
- Striving for specific LED readings or tones

Price: \$350  
Order # T8825  
Adapter cable price: \$20  
Order # T8820M

### Heart Rate Variability Biofeedback



**emWave2 Personal Stress Reliever:** emWave technology by Heartmath is an innovative approach to improving wellness and facilitating personal growth based on learning to change your heart rhythm pattern to create coherence. Used just a few minutes a day, this simple-to-use technology helps you transform feelings of anger, anxiety or frustration into more peace, ease and clarity.

Price: \$199  
Order#: emWave2



**Inner Balance:** The Inner Balance Sensor for iOS is an innovative approach to improving wellness through training, education and self-monitoring. With the Inner Balance app you enter how you are feeling in the moment and then run a session focusing on a positive moment. Research has shown that focusing on positive emotions like appreciation can greatly reduce the effects of stress and deliver you to a more calm peaceful state.

Price: \$129  
Order#: Inner Balance

Simply plug the HeartMath Sensor for iOS® onto your iPhone®, iPad® or iPod® Touch and connect the ear piece to your ear. Tap in the space below the breath pacer to start your session. Follow the pacer and the prompts below to help you moderate your breathing, activate a positive emotion and restore your Inner Balance.

## Stand Alone Instruments

### Temperature Biofeedback Trainers

**Mindfield eTemp** The eSense Temperature is a small sensor for measuring your skin temperature with the microphone input of your smartphone or tablet (Android or Apple iOS.) A free app by Mindfield® is included to display your skin temperature in different ways for biofeedback training. You can export the measured values via e-mail.



Price: \$99  
Order# ETemp

### **GSR/TEMP 2X**

The GSR Temp 2X® also allows you to do “hand warming” biofeedback. In addition, the GSR Temp 2X® monitor, the system includes a temperature sensor, remote electrodes that attach to the hands or feet, dual sensitivity meter, earphone, instruction manual and a cassette with a short relaxation program.



Price: \$140  
Order# T2120M

**SC911:** Home temperature trainer. High Contrast LCD display. Displays temperature to the nearest 0.1. Thermistor wire length - 10 ft. Monitors both body temperature and room temperature. Easy to use IN/OUT Switch. Temperature updates every 2 seconds, Accuracy +/- 1.80 F, Displays Fahrenheit or Celsius Temperature. Electronic Quartz Accuracy. Temperature range 58 F to +158 F. Features - Digital Clock, Desktop stand that tilts and thus allows for easy reading, Battery Life approximately one year. Battery included. Unit size 2 3/4 inches long, 2 1/8 inches high, 5/8 inches wide. Fits easily in shirt pocket or purse.




Price: \$19.95  
Order# SC911

For special offers and discounts, please call us at 877-669-6463

## Stand Alone Instruments


### Temperature Biofeedback Trainers

**Stress Dot Training Card:** Minimum order of 25. 10 dots per card. Match color of stress dot with carry card.  
 Price: \$1.81 Order # SDTC-0928



#### STRESS CARDS 500

Order# SC-0931  
 Price: \$1.25 EACH  
 Minimum imprinted order is 250  
 Setup charge per side/color Price: \$30  
 Logo reduction, if necessary Price: \$8




STRESS CARDS 250  
 Order# SC-093

Price: \$1.49 EACH

#### BIO DOTS (100)

Small liquid crystals circles with adhesive back  
 Price: \$21  
 Order # BD-0936



#### FINGER THERMOMETERS

Lowest cost, simply tape to end of finger. Uses non-toxic oil based liquid. No mercury. Excellent promotional piece.

Price \$1.10 EACH  
 Order# FT-0923



### Skin Conductance Biofeedback Trainers

**The GSR 2™** precisely monitors your stress levels by translating tiny tension-related changes in skin pores into a rising or falling tone.

Price: \$60  
 Order #T2001M



## Stand Alone Instruments

**Mindfield eSense** The Mindfield eSense Skin Response is a small sensor to measure your skin conductance using two electrodes and the microphone input of your smartphone or tablet (Android or Apple iOS.)



Price: \$99  
Order# ESense

### HEM 705CP SEMI AUTOMATIC BLOOD PRESSURE UNIT

Keep accurate records of your blood pressure readings with the Blood Pressure Monitor and Print Out. Simply wrap the cuff around your arm and press START. In seconds, your blood pressure and pulse are displayed on the large digital panel. You can then print your readings in either numerical or bar graph form for easy record keeping.



Price: \$190  
Order# HEM-705CP  
Adapter  
Price: \$30  
Order# HEM-0917

**Lumiscope Digital Blood Pressure:** Alternately displays blood pressure and pulse readings. Features: Touch pad control, Memory recall, Auto-off. Includes cuff (fits arms of 8.6"-12.6" in circumference). Price: \$ 45 Order# LAD-1060  
Price: \$55 Order# LADL-1060




### Respiration Feedback Devices

#### VOLDYNE INSPIROMETER

Use to practice diaphragmatic breathing. Measures volume of air inhaled.

Price: \$18 each  
Order# VI-0941  
Case of 12 \$190  
Order# VI-0942



Fully-automatic inflation upper arm blood pressure monitor  
large digital display  
Order # LAD1130  
Price \$95.00



RESP@RATE interactive therapeutic breathing device. Breathing correctly reduces unhealthy stress responses.  
Order # RES-001  
Price \$299

For special offers and discounts, please call us at 877-669-6463

## Microcurrent Pain Management

---



### The CES Ultra

Get Relief from [Anxiety](#), [Depression](#) and [Insomnia](#) — Without Medication

**The CES Ultra** is a medical device for the treatment of stress and stress-related disorders. It employs gentle electrical stimulation to the head (cranial electrical stimulation) to help normalize brain functioning.

The CES Ultra is registered with and regulated by the FDA (Food and Drug Administration).

CES units have been used internationally since the 1950s and are well-researched for their safety and effectiveness.

Often effective where drugs have failed

Portable and easy to use

Pediatric to geriatric applications

More cost-effective than medications

Non habit-forming or addictive

The CES Ultra is the only CES unit currently on the market with the proven technology and configuration employed in more than 100 scientific and research studies, including the Meta-Analysis conducted by the Harvard School of Public Health.

Price: \$349  
Order# CES-U

---

For special offers and discounts, please call us at 877-669-6463

## Software—Applications

### BIOFEEDBACK TUTOR

Developed by Biosource Software, provides an up-to-date and comprehensive review of biofeedback for BCIA certification preparation and professional education.

The 18 tutorials summarize key findings from critical articles on the BCIA reading list and major references in anatomy and physiology, biofeedback, chronic pain, clinical efficacy, heart rate variability (HRV), instrumentation, neurofeedback, physiological psychology, professional ethics, psychopharmacology, relaxation procedures, stress management, and surface electromyography.

Price: \$75                      Order# Bio-T0727

### BIOFEEDBACK 100

Developed by Biosource Source, this web-based service provides an up-to-date and comprehensive review for the BCIA Certification Exam in Biofeedback. It allows you to take 20-question exams and receive immediate feedback. The question pool exceeds 50 questions and is revised quarterly.


Price: \$50                      Order# Bio-100-0700



### NEUROFEEDBACK 100

Developed by Biosource Software, this web-based service provides an up-to-date and comprehensive review for the BCIA certification exam in EEG Biofeedback. It allows you to take 20-question exams and receive immediate feedback. The question pool exceeds 400 questions and is revised quarterly.


Price: \$50                      Order # Neuro-100-0710



### DynaMAP™ Suite

The DynaMap Suite is focused on SEMG with a wide variety of different applications for muscle re-education and relaxation. It also incorporates other modalities, including Goniometer and Force, in some of the channel sets and screens. The Suite is designed to be used on all hardware systems which support BioGraph Infiniti software with channels sets to maximize the features of each.

Price \$295                      Order# DYN-010



### BioGraph®-X Infiniti

BioGraph-X Infiniti CD 1 (83 animations) and CD 2 (55 animations) are the first in a series of animation CD's to be offered by Thought Technology to enhance your training sessions. The CD's are compiled of a variety of animations from nature, to fractals and games.

Price: \$99                      Order# BIOGX-001



## Software-Applications

### CARDIOPRO Infiti

**CardioPro™** CardioPro Infiti is a specialized off-line analysis module which brings advanced heart rate variability (HRV) analysis capabilities to the BioGraph Infiti platform. It complements any physiological biofeedback suite (such as the Physiology suite) by providing sophisticated data analysis methods that were only previously available in very specialized software packages such as Thought Technology's CardioPro software.



Price: \$600  
Order# SA7590



Price: \$295  
Order# SA7970

Physiology Suite - SA7970

#### Quick Overview

The Physiology Suite includes all the tools you need for multi-modality physiological assessments and over 80 display screens for training with standard physiological sensors including heart rate with EKG or blood volume pulse (BVP), electromyography (EMG), skin conductance, temperature and respiration. Using the Physiology suite, you can train self-regulation with physiological biofeedback protocols that include relaxation with EMG, temperature or skin conductance, heart rate variability with respiration and EKG or BVP and abdominal/thoracic breathing. The suite also features a full-length relaxation exercise script with a recorded audio track guiding the client through deep relaxation and progressive recovery



Price: \$295  
Order# SA7950


EEG Suite - SA7950

#### Quick Overview

The EEG Suite includes new coherence training screens with dual-threshold range training. It offers over 80 display screens for many standard EEG protocols, including an outstanding alpha-theta training screen and a number of template-style screens for training with 3, 6 or 10 user-defined EEG bands. The standard protocols can be used for training alpha amplitude (range), alpha peak frequency, alpha theta, beta amplitude, SMR amplitude, theta beta, theta SMR and wide band inhibit.

### T.O.V.A. –Visual & Audio

Screening tool for ADD. 22.5 min. computerized assessment used in conjunction with classroom behavior ratings. Includes 5 tests.



Visual/Audio  
Price: \$375  
Order#  
TOVA0711  
Additional tests,  
each  
Price: \$15  
Order#  
#TOVA0713




## Supplies & Electrodes

### General Supplies

#### TD-30 PREP PADS

Easiest and quickest way to prepare the skin for EMG and EEG hook-ups. Lower resistance measures. Alcohol pad impregnated with pumice Box of 100  
 Price: \$15.50  
 Order# TD-30



#### SENSOR JEL TUBE

Price \$4  
 Order#



#### NUPREP SKIN PREPPING PASTE

4 oz. Bottle  
 Price: \$30 (3 Pack)  
 Order# N-1025



#### TEN/20 CONDUCTIVE PASTE

Adhesive EEG electrode paste  
 4 oz.  
 Price: \$27 (3 Pack)  
 Order# 1020



#### Micropore™ TAPE

For Temp lead Attachment, 1 roll  
 Price: \$ 2.50  
 Order# TP-1007

12 Rolls  
 Price: \$24  
 Order# TP-1008



#### SENSOR ELECTRODE CREAM,

5 oz. Bottle  
 Price: \$4  
 Order# 1030



#### NEUROTRODE NET3

Disposable Pre gelled Hard Plastic Disc Electrode

60 electrodes per bag / 6 electrodes per foil pouch / 10 foil pouches per bag. Price: \$34.00 Order# NET3



RINGTRODE 2RT3 Disposable Pre gelled Hard Plastic Disc Electrode.

150 electrodes per bag  
 Price: \$45 Order# 2RT3



A 10/20 Universal EEG fixture (part # 10/20UF) Suggested retail \$ 75.00. For single electrodes placement, a Head band (part # HBEG 7-2 ) Suggested retail \$ 10.00 (see onset on main picture). Both fixtures allow the sensors to be tightly secured on the scalp and mastoid. 2- A cable (part # CEG 7- 06) allows connection between the sensor and the equipment, suggested retail \$ 50.00. It is compatible with most EEG Equipment.

## DISPOSABLE SENSORS

### PREGELLED DISPOSABLE ELECTRODES

Box of 50  
 Price: \$40  
 Order# DBF3D77



### QUICK STICK TYPE PERFORATED DISPOSABLE SENSORS

Package of 100 in strips of 3  
 Price: \$80  
 Order# 3BF3



### White Cloth Moveable Pre-gelled Disposable EMG Sensors

Package of 150  
 Price: \$45  
 Order #3S3GN



**Triode**, electrode with standard 2cm spacing of silver silver chloride electrodes, backed with nickel plated brass snaps to prevent corrosion when connected to pre-amplifiers for extended periods.

Price: \$85.00  
 Order # T-3402M



**UniGel** electrodes, pre gelled single electrodes, or sensitive placements on dry skin. Package of 100  
 Price: \$35  
 Order # T3425



**Ringrode** Pre-Gelled Hard Plastic Disc  
 Electrode: 50 blue ring electrode snaps per bag / 100 gel discs per bag  
 Price: \$34  
 Order # 2RT3



**UniGel** Pregelled Hard Plastic Disc Electrode: Quantity 60 (6 per foil pouch, 10 pouches)  
 Price: \$42  
 Order # NET

## Supplies & Electrodes

---

### Continence Supplies

#### **PERRY SENSORS (VAGINAL)**

Price: \$35  
Order# P6330



TTL Vaginal Sensor 38 in. 97 cm  
For Thought Technology EMG instruments  
Price: \$50  
Order #T6050

(RECTAL)  
Price: \$59.95  
Order# P6340



TTL Rectal Sensor 38 in, 97 cm  
For Thought Technology EMG Instruments  
Price: \$50  
Order #T6051



For special offers and discounts, please call us at 877-669-6463

---

## Supplies & Electrodes

### THOUGHT TECHNOLOGY SUPPLIES

#### MYOSCAN-PRO EMG SENSOR

A pre-amplified surface electromyography sensor used with the ProComp Infiniti channels C – H for RMS sEMG. Features a range switch in the sensor head to change filter settings between 0 – 400  $\mu$ V Narrow filter, 0 – 400  $\mu$ V Wide filter, 0 – 1600  $\mu$ V Wide filter. Compatible with Triode electrodes or extender cables for wider placement of electrodes.

Price: \$395 Order# T9401M



#### MYOSCAN EMG SENSOR

A pre-amplified surface electromyography sensor used with the ProComp Infiniti channels A – B for RAW sEMG. Compatible with Triode electrodes or extender cables for wider placement of electrodes. Range of 0 – 2000  $\mu$ V.

Price: \$295 Order# SA9503P



#### TEMPERATURE SENSOR

A skin temperature sensor, measuring between 10°C – 45°C ( 50°F - 115°F ).

Price: \$ 210 Order# SA9310M



#### RESPIRATION SENSOR

An easy fitting high durability latex rubber band fixed with Velcro respiration belt for monitoring respiration rate from a client. Can be worn either thoracically or abdominally, over clothing.


Price: \$ 295 Order# SA9311M



#### EMG MyoScan-Z™ Sensor - T9503Z

The new MyoScan-Z EMG sensor joins our line-up of physiological sensors. Designed for a sampling rate 2048 samples per second, this new EMG sensor offers a unique impedance check mode, which can be triggered from BioGraph Infiniti software version 5.0. It measures raw EMG signals with a range from 0 up to 2000 micro-Volts.

Price: \$395 Order# T9503Z



#### SENSOR EEG-Z3

With impedance checking: regular EEG mode, low frequency mode (0.01 Hz – 1 kHz) for evoked potentials and low frequency applications and DC mode (0 Hz – 1 kHz) for slow cortical potentials.

Price: \$495 Order # T7680M

#### HR/BVP SENSOR

A pulse detection sensor housed in a small finger worn package, to measure pulse rate of client.



Price: \$295 Order# SA9308M

#### SKIN CONDUCTANCE SENSOR

A sensor to measure the conductance across the skin, normally connected to the fingers or toes. Supplied with two finger bands.

Price: \$295 Order# SA9309M



#### EEG-Z SENSOR with Impedance checking

Price: \$395 Order # T9305Z



#### EKG SENSOR

Price: \$295 Order # T9306M



#### SENSOR EEG-Flex-Pro

Price: \$25 Order # T9305M



#### EKG SENSOR

With wrist straps Price: \$325 Order # T9307M

#### VOLTAGE ISOLATOR SYSTEM (FOR CONNECTING 3RD PARTY DEVICES) Price: \$600 Order # T9405

#### TT-pIR HEadGear™ System - T2600

The TT-pIR HEadGear features two highly accurate passive infrared (pIR) sensors in a wide-range, dual configuration providing stable and reliable measurements. The lightweight, ergonomic design provides maximum comfort and an open concept that greatly reduces ambient temperature dependence and stabilization time.

Comes equipped with TT-pIR Mini-Suite software for monitoring, training and follow-up.

Price: \$ 295 Order# T2600

## Audio CD's, & DVD Video

CHANGE THE CHANNEL ON PAIN, By Emmett Miller,  
MD

Price: \$17.95                      Order# SC5046    CD

**DEEP MUSCLE RELAXATION**

Price: \$16.95    Order #CL-101HCD    CD

**Self Directed Relaxation (Autogenic Relaxation)**

Price \$16.95    Order # 102HCD    CD

**FREEING YOURSELF FROM FEAR**

CD \$39.95            Order # SC5338                      3 CDs

**HEADACHE RELIEF**

By Emmett Miller, MD

Price: \$16.95    Order # SC5012                      CD

**DOWN WITH HIGH BLOOD PRESSURE**

Price: \$16.95    Order# SC5003                      CD

**FOCUSED MIND STATE**

Maximizing Your Potential Through the Power of Con-  
centration (6 CD)

Price: \$69.95

Order # FMS-0207CD                      CD

**HEALING JOURNEY**

By Emmett Miller, MD

Price: \$16.95    Order # SC5016                      CD

**HEALTH AND WELLNESS**

By Emmett Miller, MD

Price: \$16.95    Order # SC5018                      CD

Gentle Mood Music (solo guitar) for Background for Relaxa-  
tion

Order# HHCD                      Price: \$16.95

Relax... Let Go ... Relax (Scan body for tension – then dis-  
charge through breathing and releasing)

Audio CD            Order# HRLXCD Price: \$16.95

**StressRelease**

Audio CD            Order# HSRLCD                      Price: \$16.95

**Eating**

Audio CD            Order# HECD                      Price: \$16.95

**Smoke No More**

Price: \$25.50    Order # SC5221 2 CDs

**LETTING GO OF STRESS** - World's Most Popular  
Stress Tape Using Guided Imagery Musical Score by Steve  
Halpern

Price: \$16.95                      Order # SC5023    CD

**Gateway to Inner Peace**

For Deep Relaxation and Self-Healing

by Sarena S. Morello, M.S., LMHC (Licensed Mental  
Health Counselor The program comes with: Audio CD,  
Biodots for biofeedback ,5 page instruction booklet

Price: \$18            Order # GIP001

**OPEN FOCUS**, By Dr. Les Fehmi

6 Tapes—Intro, Head/hands, Long & Short Form, Dissolv-  
ing Pain, Independent Practice

Price: \$79.95    Order# OF0231CD                      CDs

**RELAX**

By Dennis Weaver. VHS Tape.

Price: \$29.95    Order# R0214

**RELAXATION & INSPIRATION**

Price: \$16.95    Order # SC5102                      CD

**SELF-GUIDED RELAXATION** By JOEL LEVY

Price: \$17.95    Order# SGR0221

**TEN MINUTE STRESS MANAGER**

By Emmett Miller, A favorite of the U.S. Olympic Team!

Price: \$16.95    Order # SC5053 CD

**STRESS CONTROL FOR KIDS** - 6 Audio CDs Price:

\$30    Order# 200HS

## Audio CD's, & DVD Video

---

Stressfitness Vol. I, Emmett Miller

Order #SC5021 (1 CD) Price \$16.95  
Stress Fitness Vol. II

Order #SC5022 (1 CD) Price: \$16.95

TIME TO RELAX, By Dr. Lowenstein

Price: \$16.95 Order # CL115HC CD

SLEEPEAZY Fall Asleep, Stay Asleep, Easy™

by Sarena S. Morello, M.S., LMHC

program includes: Videotape, Cassette tape, 36-page  
booklet, Biofeedback finger temperature trainer

Price: \$39.95 Order # SE021

Stress Inoculation Video, by Matthew McKay, Ph.D.

Price: \$49.95 Order# 318X

Facial Paralysis: Neuromuscular Re-Education with sEMG Bio-  
feedback CD-ROM, Gabriel Sella

Price: \$40 Order# BC-63

Myofascial Pain Syndrome: Manual Trigger Point & sEMG Bio-  
feedback Therapy Methods CD-ROM

Price: \$30 Order# BC-64

Stress Release

Audio CD Order# HSRLCD Price: \$16.95

Body Image

Audio CD Order# HBICD Price: \$16.95

Warm and Heavy (Classic Autogenic Relaxation)

Audio CD Order# HWHCD Price: \$16.95

Spiritual Centering

Audio CD Order# HSPCCD Price: \$16.95

Healthy Balancing

Audio CD Order# HHBTCDD Price: \$16.95

Inner Healing

Audio CD Order# HHHCD Price: \$16.95

Personal Empowering

Audio CD Order# HPECD Price: \$16.95

Countdown to Relaxation

Audio CD Order# HCDCD Price: \$16.95

Take a Deep Breath

Audio CD Order# HTDBCDD Price: \$16.95

Mantras

Audio CD Order# HMTCD Price: \$16.95

Daydreams 1 Getaways

Audio CD Order# HDD1CD Price: \$16.95

Daydreams 2 Peaceful Places

Audio CD Order# HDD2CD Price: \$16.95

Daydreams 3 Relaxing Retreats

Audio CD Order# HDD3CD Price: \$16.95

Breaktime, Stress Break

Audio CD Order# HBTCDD Price: \$16.95

Natural Tranquilizers

Audio CD Order# HNTCD Price: \$16.95

Stress Escapes

Audio CD Order# HSECD Price: \$16.95

Healthy Choices

Audio CD Order# HHCCD Price: \$16.95

Healing Visions

Audio CD Order# HHVCD Price: \$16.95

Refreshing Journeys

Audio CD Order# HRJCD Price: \$16.95

---

## Audio CD's, & DVD

Wilderness Daydreams 1 Canoe/Rain Audio CD Order# HWD1CD Price: \$16.95	Video Library Price: \$99 each
Wilderness Daydreams 2 Island/Spring Audio CD Order# HWD2CD Price: \$16.95	Making Healthy Choices
Wilderness Daydreams 3 Campfire/Stream Audio CD Order# HWD3CD Price: \$16.95	Healthy Lifestyle HMHV1
Massage Audio CD Order# HMACD Price: \$16.95	Healthy Eating HMHV2
Yoga Audio CD Order# HYCD Price: \$16.95	Healthy Exercise HMHV3
Worry Stoppers Audio CD Order# HWSCD Price: \$16.95	Healthy Stress HMHV4
Calm Down Audio CD Order# HCALMCD Price: \$16.95	Healthy Relationships HMHV5
Becoming Light-Hearted, Use Your Sense of Humor to Manage Stress (audio CD) Order# HBLH Price: \$16..95	Healthy Change HMHV6
MUSIC ONLY – NO VOICE – Steven Halpern, Ph.D. Contemplation Audio CD Order# HCTCD Price: \$16.95	Workplace Relationships Building A Positive, Productive Environment HWRV1
Harmony Audio CD Order# HHCD Price: \$16.95	Negotiating Conflict HWRV2
<b>RELAXATION VIDEO AND DVD</b>	Preventing Sexual Harassment HWRV3
The Healing Power of Nature Video Tape Order# HPNV Price: \$29.95	Addressing Diversity HWRV4
Natural States Video Tape Order# HNSV Price: \$29.95	Working With Difficult People HWRV5
America the Beautiful Video Tape Order# HABV Price: \$29.95	Meeting the Challenge of Change HWRV6
Earth Dance Video Tape Order# HEDV Price: \$29.95	Managing Job Stress
The Art of Nature Video Tape Order# HART Price: \$29.95	Handling Workplace Pressure HJSV1
Take the Boundary Waters Home With You DVD Order# HBWDVD Price: \$19.95	Clarifying Roles and Expectations HJSV2
The Seasons of the North Shore DVD Order# HSEASONS Price: \$19.95	Controlling the Workload HJSV3
Alaska: Mystic Ice	Managing People Pressures HJSV4
	Surviving the Changing Workplace HJSV5
	Balancing Work and Home HJSV6
	Manage It
	Stress Traps HMIV1
	Stress Overload HMIV2
	Interpersonal Conflict HMIV3
	Addictive Patterns HMIV4
	Job Stress HMIV5
	Survival Skills HMIV6
	Success with Stress
	Take Charge of Your Life HSSV1
	Get Ready for Change HSSV2
	Build Healthy Relationships HSSV3
	Manage a Fast-Paced Life HSSV4
	Replace Self-Defeating Behaviors HSSV5
	Develop New Stress Skills HSSV6
	Humor Works, Izzy Gesell A Fresh and Creative Approach to Stress Management (A 3-Session Video Course) HIHV

## Books

---

### **MINDFITNESS TRAINING:NEUROFEEDBACK AND THE PROCESS**

By R. Adam Crane, Richard Soutar  
Price: \$21.95      Order# 0100

### **BIOFEEDBACK A PRACTITIONER'S GUIDE**

Schwartz & Associates  
Price: \$65      Order# 0103

### **BIOFEEDBACK & SOMATICS**

By Eleanor Criswell  
Price: \$14.95      Order# 0104

Getting Rid Or Ritalin: How Neurofeedback can successfully Treat Attention Deficit Disorder Without Drugs

By Robert W. Hill  
Price: \$19.95      Order# 1571742549

### **HUMAN BRAIN COLORING BOOK**

By DRS. Diamond, Scheibel, Elson  
Price: \$16      Order# 0111

### **INTRODUCTION TO QUANTITATIVE EEG AND NEUROFEEDBACK**

By James R. Evans, Andrew Abarbanel ed.  
Price: \$108.95      Order# 0112

### **The Open Focus Handbook**

By Drs. Fritz and Fehmi  
Price: \$29.95      Order# 0124

### **Handbook of Neurofeedback**

By James R. Evans  
Price: \$99.95      Order#0124B

### **STRESS PASSAGES**, By Dr. Mason

Price: \$11.95      Order# 0118

### **A SYMPHONY IN THE BRAIN**, By Jim Robbins

Price: \$13.50      Order# 0122

**Relaxation and Stress Reduction Workbook**, by Martha Davis, Ph.D., Elizabeth Robbins Eshelman, MSW, and Matthew McKay, Ph.D. Price: \$19.95      Order #2140

Getting Started in Neurofeedback, John Demos  
Price: \$35      Order# GSN-001

**WISDOM AT WORK**, By Dr. Joel Levy and Michelle Levy  
Price: \$13.95      Order# 1-57324-180-6

**The Neurofeedback Book**, Michael and Lynda Thompson  
Price: \$99      Order # BC-76

**What Stress Can Do**, Why you must recognize all of the damage stress is causing in your life and do something about it, By Harry L. Campbell

Price \$14.95      Order# HLC-001



## BCIA Certification Training

### Professional Peripheral Biofeedback BCIA Program

This innovative seminar is designed to teach clinicians biofeedback fundamentals and cutting-edge applications. This program covers all 48 hours of the BCIA Didactic Biofeedback Education requirement through powerful demonstrations of biofeedback and adjunctive procedures, imaginative hands-on lab exercises with multi-channel computerized equipment, professional multimedia presentations, and comprehensive review materials. This program integrates personal hands-on training with equipment with didactic presentation.

Participants will learn:

The foundations of biofeedback

Physiological basis and recording of surface electromyography (sEMG), temperature, electrodermal activity (EDA), respiration patterns (strain gauges and pCO<sub>2</sub>), electroencephalography (EEG), peripheral blood flow (BVP), heart rate and respiratory sinus arrhythmia (RSA), and economic biofeedback approaches (including incentive spirometers and scales)

How to identifying artifacts and determining appropriate recording

How to record and interpret data

How to perform a psychophysiological profile and pre and post data comparison

The use of adjunctive techniques such as progressive relaxation training, autogenic training, imagery, guided breathing, desensitization, mindfulness training, quieting reflex.

The role of home biofeedback devices such as Freezeframer, Wild Divine, and GSR2

Overview of clinical and educational applications

Price: \$1,195 #BCIA-B

### Professional EEG Neurofeedback Certification Program

This competency based curriculum includes practical, easily learnable strategies which may be used with virtually all of the most common EEG Neurofeedback applications. Designed to give the practitioner the knowledge and tools needed in order to rapidly achieve clinical and economic success in EEG Neurofeedback, this program emphasizes those neurofeedback adjunctive applications for which there is the greatest demand, including but not limited to: ADD / ADHD / Conduct Disorder & other Learning disabilities, Autism Disorders, Substance Abuse, eating disorders and addictive behaviors, PTSD, Multiple Personality Disorder, Depression, Behavioral and Stress Related Disorders, Anxiety / Panic Attacks, Pain Management, Chronic Pain, headaches (Migraine and Tension), RSD, Fibromyalgia, Phantom Limb etc., Adult and childhood sleep disorders, Traumatic Brain Injury, Stroke Rehabilitation, Dementia (non-Alzheimer's), Epilepsy, Chronic Fatigue Syndrome and MindFitness® Life / Performance Enhancement. Includes demonstrations and hands-on training on the most popular, versatile, cost effective, reliable and well supported EEG neurofeedback systems made as well as introductions to neurocardiology (Heart Rate Variability), Temperature, EMG and other complementary modalities, lively dialogues, case histories and the information you need in order to select the right equipment for your circumstances.

Many practitioners have told us that they had attended other EEG neurofeedback seminars but this is the one that put it all together for them and showed them how to provide quality neurofeedback services while clarifying many theories and hypotheses associated with EEG Neurofeedback.

Previous experience in biofeedback, EEG neurofeedback or computers is not required.

Supervised evening sessions provide additional opportunity for hands on equipment training.

Topics in the course include:

Orientation to EEG Biofeedback

Basic Neurophysiology and Neuroanatomy

Treatment Planning

Instrumentation & electronics

Research

Psychopharmacology Considerations

Professional Conduct

Price: \$995 #BCIA-N

For special offers and discounts, please call us at 877-669-6463

## Computer Peripherals

---



A partial description of some of the hardware that we offer is:

**COMPUTERS:** Desktops, Towers, Notebooks

**MONITORS:** Standard, flat screen, 17", 19", 21".

**PRINTERS:** Laser, Inkjet

Brand names include but are not limited to: Toshiba, Dell, Viewsonic, Epson, Hewlett Packard



Please see our training brochure for information on our BCIA Certification Biofeedback, EEG/Neurofeedback

**SUPPORT:** We can answer your clinical, software, and technical questions. We can assist you in choosing the right components. We help you install and run hardware and software quickly and easily.

**Consulting, Training Onsite or By Phone:** \$95 per hour

**Onsite Consulting and Training:** \$660 - \$1,700 per day (Price depends on need)



### Terms and Conditions

Freight, Shipping, and Handling added to listed prices. Minimum charge is \$8.00 per shipment. Visa, Master Card, and American Express are accepted. All prices F.O.B. All Prices are in U.S. Dollars. Leasing is available.

For International Payment Terms: Foreign payment may be made in advance by Visa, Master Card, Wire transfer of funds to Biofeedback Resources' bank account or by confirmed irrevocable letter of credit. (All banking charges and fees must be applied to the buyer's account). The validity period of letters of credit should be at least 90 days. All funds should be in US dollars.

Orders paid by direct transfer of funds or by credit card can be shipped via air express, standard air freight, or via air mail (some size restrictions may apply to air mail shipments). Orders placed with a letter of credit will be shipped by standard air freight. Please contact us for a pro forma quotation to include shipping and any additional costs before placing an order.

Return Policy/Refund Policy: Products may be returned within 15 days of receipt of order for a refund of the purchase price less a 15% restocking fee and shipping and handling charges. Defective products may be returned within 10 days of receipt for a refund or replacement.

Price list effective June 1, 2015 Prices subject to change without notice.



109 Croton Ave.  
2nd Floor, Suite #5  
Ossining, NY 10562